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Interferencias: Heparina

Referencias

1. JJM van Dongen et al. Standardized RT-PCR analysis of fusion Gene Transcripts from chromosome aberrations in acute leukemia for detection of minimal residual disease. *Leukemia* 1999 (13) 1901-1928.
2. Niels Pallisgaard et al: Multiplex reverse transcription-polymerase chain reaction for simultaneous screening of 29 translocation and chromosomal aberrations in acute leukemia. *Blood* 1998 (92) 574-588.
3. Marta Libura et al. FLT3 and MLL intagenic abnormalities in AML reflect a common category of genotoxic stress. *Blood* 2003 (102) 2198-2204.
4. Baxter et al. Acquired mutation of the tyrisine kinase JAK2 in human myeloproliferative disorders. *Lancet* 2005 (365) 1054-61.
5. Peter W M Johnson. Polymerase chain reaction for detection of the t(14;18) translocation in Lymphomas. *Methods in Molecular Mediicine. Molecular diagnosis of Cancer*, FE Cotter Humana Press Inc. Totowa NJ.
6. Molot R J et al. Antigen expression and polymerase chain reaction amplification of mantle cell lymphomas. *Blood* 1994 (83) 1626-31.
7. Yee H T et al. Molecular characterization of the t(2;5)(p23; q35) translocation in anaplastic large cell lymphoma (Ki-1) and Hodgkin's disease. *Blood* 1996 (3) 1081-88.
8. Barr F G et al. Genetic heterogeneity in the alveolar rhabdomyosarcoma subset without typical gene fusions. *Cancer Research* 2002 (62) 4704-10.
9. Zucman J et al. Combinatorial generation of variable fusion proteins in the Ewing family of tumors. *EMBO J* 1993 (12) 4481-87.
10. Hiroyuki N et al. Detection of tyrosine hydroxilase mRNA and minimal neuroblastoma cells by the reverse transcription-polimerase chain reaction. *Eur J Cancer* 1991 (27) 762-765.
11. Gilbert J et al. Determination of N-myc gene amplification in neuroblastoma by differential polymerase chain reaction. *Mol Cell Probes* 1993 ; 7(3) 227-34.
12. Schaumacher J A et al. Detection of the c-kit D816V mutation in systemic mastocytosis bya allele-specific PCR. *J Clin Pathol* 2008 (61) 109-114.
13. Sotlar K et al. Detection of c-kit mutation Asp 816 to Val in microdissected bone marrow infiltrates in a case of systemic mastocytosis associated with chronic myelomonocytic leukaemia. *J Clin Pathol* 2000 (53) 188-93.
14. Gorello P et al. Quatitative assessment of minimal residual disease in acute myeloid leukemia carrying nucleophosmin (NPM1) mutations. *Leukemia* 2006 (20) 1103-1108.
15. Brunangelo F et al. Cytoplasmic nucleophosmin in acute myelogenous leukemia with normal karyotype. *New England J of Medecine* 2005 (352) 254-266.
16. Gedeon L et al. FIP1L1-PDGFRA molecular analysis in the differential diagnosis of eosinophilia. *Blood Disord* 2009 (9).
17. Wei X et al. MPL W515L mutation in chinese patients with myeloproliferative diseases. *Leukemia & Lymphoma* 2008 (49) 955-958.
18. Brandford S et al. High frequency of point mutations clustered within the adenosine triphosphate-binding region of BCR-ABL in patients with chronic myeloid leukemia or Ph- positive acute lymphoblastic leukemia who develop imatinib (STI571) resistance. *Blood* 2002 (99) 3472-3475.
19. Hughes T et al. Monitoring CML patients responding to treatment with tyrosine kinase inhibitors: review and recommendations for harmonizing current methodology for detecting BCR-ABL transcripts and kinase domain mutations and for expressing results. *Blood* 2011 8108) 28-36.
20. Martine P et al. A multiplex real time PCR assay for the detection of gene fusion observed in solid tumors. *Laboratory Investigation* 2001 881) 905-12.
21. Gabert J et al. Standardization and quality control studies of real time quantitative reverse transcriptase polymerase chaun reaction of fusion gene transcripts for residual disease detection in leukemia- A Europe Against Cancer Program. *Leukemia* 2003 (17) 2318-2357.
22. Beillard E et al. Evaluation of candidate control genes for diagnosis and residual disease detection in leukemic patients using real time quantitative reverse-transcriptase polymerase chain reaction (RQ-PCR)- A Europe Against Cancer Program. *Leukemia* 2003 (17) 2474-2486.
23. Klampfl Tet al. Somatic Mutation of Calreticulin in Myeloproliferative Neoplasms. *N Engl J Med* (2013) 369; 25.