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Referencias

1. College of American Pathologist (1989). Clinical Laboratory Handbook for patient preparation and specimen handling. Fascicle V.
2. Mayo Laboratories. 17-Hydroxyprogesterone, Serum. Mayo Clinic.
3. Siemens, Coat a Count 17  $\alpha$ -OH Progesterone RIA kit. Instruction Manual, última versión.
4. Westgard J. Biologic Variation Database. <http://www.westgard.com/biodatabase1.htm>.
5. L. Kyriakopoulou, M. Yazdanpanah, D.A. Colantonio, M.K. Chan, C.H. Daly, K. Adeli. A sensitive and rapid mass spectrometric method for the simultaneous measurement of eight steroid hormones and CALIPER pediatric reference intervals. Clinical Biochemistry 46 (2013) 642-651
6. Therese Koala, Diane Schmiederera, Hai Pham-Tuana, Cornelia Röhring, Manfred Rauh Standardized LC-MS/MS based steroid hormone profile-analysis. Journal of Steroid Biochemistry & Molecular Biology 129 (2012) 129- 138.
7. B.K. Matuszewski, M.L. Constanzer, C.M. Chavez-Eng, Anal Chem, 75 (2003) 3019-3030.
8. FDA: Guidance for industry-bioanalytical method validation. [www.fda.gov/downloads/Drugs/.../Guidances/ucm070107.pdf](http://www.fda.gov/downloads/Drugs/.../Guidances/ucm070107.pdf)
9. Clinical Laboratory Standard Institute (CLSI) Guideline. C50-A Vol27N°3; Mass Spectrometry in the Clinical Laboratory: general principles and guidance.
10. Peña Adolfo. (2010). Pruebas de función de la corteza adrenal. Disponible en <http://adolfoneda.com/pruebas-de-funcion-de-la-corteza-adrenal/>
11. Clin Chem Lab Med 2020; 58(5): 741-752. James M. Hawley\*, Joanne E. Adaway, Laura J. Owen and Brian G. Keevil. Development of a total serum testosterone, androstenedione, 17-hydroxyprogesterone, 11 $\beta$ -hydroxyandrostenedione and 11-ketotestosterone LC-MS/MS assay and its application to evaluate pre-analytical sample stability
12. Clin Chem Lab Med 2016; Sophie Hepburn\*, Michael J.P. Wright, Conchita Boyder, Renee C. Sahertian, Ben Lu, Rui Zhang, Chris P. White and Andrea R. Horvath. Sex steroid hormone stability in serum tubes with and without separator gels