## MONITOREO TERAPÉUTICO DE EVEROLIMUS (Certican®) POR LC-MS/MS

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**Referencias**

- :
  1. A rapid HPLC-MS/MS method for the simultaneous quantification of cyclosporine A, tacrolimus, sirolimus and everolimus in human blood samples. Christoph Seger, Karin Tentschert, Wolfgang Stoggl, Andrea Griesmacher, Steven L Ramsay. (2009) Vol.4 N°.4, Nature protocols.
  2. Robust, High-Throughput LC-MS/MS Method for Therapeutic Drug Monitoring of Cyclosporine, Tacrolimus, Everolimus, and Sirolimus in Whole Blood. Remco A. Koster, BSc, Eli C. F. Dijkers, MSc, PharmD, and Donald R. A. Uges, PharmD, PhD. Ther Drug Monit \_ Volume 31, Number 1, February 2009.
  3. High-throughput analysis of everolimus (RAD001) and cyclosporin A (CsA) in whole blood by liquid chromatography/mass spectrometry using a semi-automated 96-well solid-phase extraction system. Louis M. McMahon\*, Suyi Luo, Michael Hayes and Francis L. S. Tse. Rapid Commun. Mass Spectrom. 14, 1965-1971 (2000).
  4. Simultaneous determination of four immunosuppressants by means of high speed and robust on-line solid phase extraction-high performance liquid chromatography-tandem mass spectrometry. Therese Koal, Michael Deters, Bruno Casetta, Volkhard Kaever. Journal of Chromatography B, 805 (2004): 215-222.
  5. Validation of an LC-MS/MS method to determine five immunosuppressants with deuterated internal standards including MPA. Armin Buchwald, Karl Winkler and Thomas Epting. Buchwald et al. BMC Clinical Pharmacology 2012, 12:2 <http://www.biomedcentral.com/1472-6904/12/2>
  6. Simultaneous Determination of Cyclosporine A, Tacrolimus, Sirolimus, and Everolimus in Whole-Blood Samples by LC-MS/MS. Mustafa Karapirli, Murat Kizilgun, Ozgur Yesilyurt, Husamettin Gul, Zeki Ilker Kunak, Emin Ozgur Akgul, Enis Macit, Tuncer Cayci, Yasemin Gulcan Kurt, Ibrahim Aydin, Hakan Yaren, Melik Seyrek, Erdinc Cakir, and Halil Yaman.

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7. Mayo Medical Laboratories. Everolimus, Blood. Mayo Clinic.
8. Stability of sirolimus and everolimus measured by immunoassay techniques in whole blood samples from kidney transplant patients. D. Capone et al. International Journal of Immunopathology and Pharmacology. (2008) Vol. 21, N° 2, 297-307.
9. Clinical Laboratory Standard Institute (CLSI) Guideline. C50-A Vol27N°3; Mass Spectrometry in the Clinical Laboratory: general principles and guidance.
10. FDA: Guidance for industry-bioanalytical method validation. [www.fda.gov/downloads/Drugs/.../Guidances/ucm070107.pdf](http://www.fda.gov/downloads/Drugs/.../Guidances/ucm070107.pdf)
11. B.K. Matuszewski, M.L. Constanzer, C.M. Chavez-Eng, Anal. Chem. 75 (2003) 3019.